

Presperse Finished Products

Enlighten Me!

(Skin Lightening Serum)

(#36-16-1)

WHAT IS IT?	A roll-on serum utilizing PRESPERSE'S VITAMIN C DERIVATIVE and BRIGHTENING BOTANICALS to lighten and wake up dull aging skin.
HOW DO YOU USE IT?	Gently massage the roll-on applicator onto dull age spots and dark under-eye circles. For entire face application, apply to the large face areas and gently pat into skin.
WHO SHOULD USE IT?	Anyone seeking a natural non-bleaching skin lightening treatment to combat signs of dull age spots and dark circles around the delicate eye area. Those who seek a gentle form of skin lightening, which is hydroquinone free. Ideal for the entire face to brighten up and improve the overall skin-tone
WHAT DOES IT DO?	Powerful combination of three skin lightening actives, brighten and awaken dull, hyper-pigmented skin. Ascorbyl PM , a powerful skin lightener and antioxidant, revitalizes the skin while effectively lightening it. Natural Mulberry Root Extract lightens problematic age spots and combats dulling skin through inhibition of both tyrosinase and melanine formation. Oligogeline , a water-soluble active, provides silicone-like slip and provides film-forming properties.
WHAT'S IN IT?	Water, Glycerine, Chondrus Crispus (carrageenan) Extract, Xanthan Gum, Butylene Glycol, Morus Alba Root Extract, Magnesium Ascorbyl Phosphate, Raffinose, Sodium Polyacrylate, PVM/MA Copolymer, Phenoxyethanol, Caprylyl, Ethylhexylglycerin, Hexylene Glycol, Mica, Titanium Dioxide, Disodium EDTA.

Trade Name	INCI Name	Supplier	% W/W
PHASE A			
Deionized Water	Water	N/A	61.85
Zilgel Oil	Glycerin (and) water (and) Sodium Polyacrylate (and) PVM/MA Copolymer	Presperse, LLC	3.00
Glycerin (99.7%)	Glycerin	Jeen International	3.00
Oligogeline	Water (and) Chondrus Crispus (Carrageenan) Extract	Presperse, LLC	2.50
Botanistat PF-64	Phenoxyethanol (and) Capryl Glycol (and) Ethylhexylglycerin (and) Hexylene Glycol	Botanigenics, Inc.	1.00
Disodium EDTA	Disodium EDTA	N/A	0.05
Phase B			
Butylene Glycol	Butylene Glycol	Ruger Inc.	0.45
Kelco CG	Xanthan Gum	CP Kelco	0.45
PHASE C			
Ascorbyl PM	Magnesium Ascorbyl Phosphate	Presperse, LLC	0.20
Butylene Glycol	Butylene Glycol	Ruger Inc.	0.50
Deionized Water	Water	N/A	25.00
PHASE D			
Oliga GGF	Raffinose	Presperse, LLC	0.20
Phycol CC-BG	Butylene Glycol (and) Water (and) Chondrus Crispus Extract	Presperse, LLC	1.00
Mulberry Root Extract	Morus Alba Root Extract	Presperse, LLC	0.50
PHASE E			
Butylene Glycol	Butylene Glycol	Ruger Inc.	0.15



P R O C E D U R E

1. Charge main kettle with water and add PHASE A ingredients in listed order, mix at medium speed until uniform (1-2 min) at 20-25°C.
2. Combine Phase B ingredients and add to PHASE A, continue mixing until a homogeneous phase free of solids is achieved (20-25°C).
3. In a separate container, combine PHASE C in order listed and mix with prop mixing at low speed for 1-2 minutes until a clear, uniform phase is achieved. Add PHASE C to PHASE AB with slow prop mixing
4. Add Phase D ingredients to PHASE ABC in order listed, with medium prop mixing at 20-25°C.
5. Combine Phase E ingredients and add to Batch, with low prop mixing at 20-25°C until a homogenous phase is achieved. Transfer to a closed container

INGREDIENT	FUNCTION
Ascorbyl PM	Vitamin C derivative. Skin lightening and free radical scavenger.
Zilgel Oil	Superb moisturizer, provides slip to formulation
Oliga GGF	Vegetable based humectant, provides skin conditioning
Phycol CC-BG	Film forming, moisturizing, skin protection.
Mulberry Root Extract	Skin lightening agent, both melanine formation inhibition and tyrosinase inhibition.
Oligogeline	Film former and moisturizer, provides silicone-like slip to formulation
PresPearls Violet	Visual Enhancement

S P E C I F I C A T I O N S

APPEARANCE (@25C)	Transparent, pearlescent, low viscosity fluid
Odor	Characteristic
pH	5-7

PRESPERSE MAKES NO WARRANTY, EXPRESSED OR IMPLIED, OTHER THAN THOSE EXPRESSLY SET FORTH ABOVE, AND ALL OTHER EXPRESSED AND IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE HEREBY EXPRESSLY DISCLAIMED BY PRESPERSE. WITHOUT LIMITING THE GENERALITY OF THE FOREGOING, PRESPERSE MAKES NO WARRANTIES, EXPRESSED OR IMPLIED, WITH RESPECT TO THE USE OR MISUSE OR REFORMULATION OF GOODS SOLD HEREUNDER ON OR WITH GOODS NOT SOLD BY PRESPERSE. PRESPERSE DOES NOT WARRANT FAVORABLE RESULTS FROM, AND ASSUMES NO LIABILITY IN CONNECTION WITH THE USE OF THE PRODUCTS SOLD HEREUNDER. THE PRODUCTS SOLD HEREUNDER ARE INTENDED FOR USE BY PERSONS HAVING TECHNICAL SKILL AND KNOWLEDGE, AT THEIR OWN DISCRETION AND RISK.

© 2009 Presperse, LLC. All rights reserved.

FORMULATION #: 36-16-1

DATE: 4/24/09